



SpinStar™ Viral NA Kit Ultra

Ultra sensitive extraction of viral
nucleic acids (DNA and RNA)

Highest sensitivity for low sample volume

The SpinStar™ Viral Nucleic Acid Kit Ultra is intended for simultaneous **ultra sensitive purification of viral nucleic acids from 200 µl** of fresh or frozen human serum and plasma, respiratory specimens in VTM and stool.

The kit is a validated one stop solution for extraction of both viral DNA and viral RNA in a single workflow and offers fast, convenient and consistent high quality and ultra sensitive results through proven silica membrane technology in spin column format. The purified DNA and/or RNA is ready to use for downstream analysis like PCR, real-time PCR or any other molecular application that requires pure nucleic acids of high quality and integrity.



Order Information

Product	Content	Reactions	Cat. No.
SpinStar™ Viral Nucleic Acid Kit Ultra 1.0	SpinStar™ Ultra columns, collection tubes, lysis buffer, wash buffer 1, wash buffer 2, carrier RNA, Proteinase K, elution buffer, user manual	100	813803

ADT Biotech

307 Block B, Phileo Damansara 1, 9 Jalan 16/11
46350 Petaling Jaya, Selangor, Malaysia

phone +603 7931 6760
fax +603 7931 5352
email info@adt-biotech.com
website www.adt-biotech.com

