



SpinStar™ VNAplus Mix

For use in combination with
SpinStar™ Viral Nucleic Acid Kit

Ultra-sensitive extraction of viral targets

The SpinStar™ VNAplus Mix allows concentration of viral DNA and RNA from 1 ml plasma or serum, prior to purification with the SpinStar™ Viral Nucleic Acid Kit. Ultracentrifugation and specialized laboratory equipment are not required.

The SpinStar™ VNAplus Mix forms complexes with nucleic acids, which can be sedimented by centrifugation to form a pellet. This pellet is then resuspended in standard volume of lysis buffer and nucleic acids are subsequently purified using the SpinStar™ Viral Nucleic Acid procedure, enabling efficient and ultra-sensitive viral nucleic acid extraction.



Order Information

Product	Content	Reactions	Cat. No.
SpinStar™ VNAplus Mix	90 ml ready-to-use Pretreatment Solution	100	812803

ADT Biotech

307 Block B, Phileo Damansara 1, 9 Jalan 16/11
46350 Petaling Jaya, Selangor, Malaysia

phone +603 7931 6760
fax +603 7931 5352
email info@adt-biotech.com
website www.adt-biotech.com

