



www.marvnier.com

ZHEJIANG MARVNIER BIOTECHNOLOGY CO.,LTD.


HIGH INTENSITY • TOP QUALITY • CERTIFIED

1.5 ML MICROCENTRIFUGE TUBE



LOCK-CAP DESIGN

Secure snap closure minimizes accidental opening



1.5 ML CAPACITY

Graduated marking for easy volume marking



PREMIUM QUALITY

Made from clear polypropylene for reliable results

DNASE/RNASE-FREE,
PYROGEN-FREE

MAX CENTRIFUGATION
FORCE (14,000–20,000 × G)

MDA-CERTIFIED

TEMPERATURE RANGE
(-80°C TO +121°C)

LEAK-PROOF CAP
DESIGN

CLEAR GRADUATION
MARKS

CONICAL BOTTOM FOR
MAXIMUM RECOVERY

LEAK-PROOF CAP
DESIGN

STANDARDIZATION OF
COMPATIBILITY

Ordering Information

Cat No.	Description	Max Volume (ml)	Packaging
MMCT-1.5	<ul style="list-style-type: none"> • DNase-, RNase-, ATP-, and pyrogen-free • With cap (attached snap-cap design) • Standard / Low-retention 	1.5	500pcs/bag, 10bags/case



| THE MARVNIER PRODUCTS